

## Leber Hereditary Optic Neuropathy (LHON) Variant Panel

**Variant List:** G3376A, G3460A, G3635A, G3697A, G3733A, G3733C, G3890A, C4171A, G10197A, T10663C, G11778A, G13042A, G13051A, T13094C, G13513A, G14459A, C14482A, C14482G, T14484C, A14495G, C14568T

### Clinical Features:

Leber hereditary optic neuropathy (LHON) is characterized by bilateral, painless subacute visual failure that develops in young adults with males being 4-5 times more likely than females to be affected.<sup>1</sup> Individuals with LHON are usually asymptomatic until developing blurred vision in the central visual field in one eye.<sup>1</sup> Similar symptoms appear in the other eye approximately 2-3 months later.<sup>1</sup> In approximately 25% of cases, both eyes are affected at onset.<sup>1</sup> After onset, the optic discs become atrophic. Significant improvements in vision are rare and most patients become legally blind.<sup>1</sup> Cardiac arrhythmias, postural tremor, peripheral neuropathy, nonspecific myopathy, and movement disorders are reported to be more common in patients with LHON than in the general population.<sup>1</sup> A multiple sclerosis-like illness has also been reported, mostly in women.<sup>1</sup>

### Genetics:

LHON is caused by variants in the mitochondrial DNA (mtDNA). Four variants (3460 G>A, 11778 G>A, 14459 G>A and 14484 T>C) account for approximately 95% of patients with LHON, while the other variants included in this panel are expected to be rare causes of LHON.<sup>1,2</sup> Each mitochondrion has multiple copies of mtDNA and there are hundreds to thousands of mitochondria per cell, dependent on the cell type. Usually mtDNA variants affect only a fraction of the mtDNA; the coexistence of normal and mutant mtDNA is called heteroplasmy. Patients with LHON generally have more than 70% mutant mtDNA in leukocytes.<sup>1</sup> Variants causing LHON have reduced penetrance; approximately 50% of males and 90% of females harboring a variant are unaffected.<sup>1</sup> Variants in mtDNA arise *de novo* or are maternally inherited. In most cases, mtDNA point variants are inherited.

### Test Methods:

Using genomic DNA, the entire mitochondrial genome is amplified by long-range PCR and sequenced for the detection of the variants on this panel using Next Generation sequencing that allows deep parallel sequencing of large numbers of mtDNA variants concurrently.<sup>14</sup> DNA sequences are assembled and compared to the published mitochondrial genome reference sequences for analysis. The presence of any variants is confirmed by conventional dideoxy sequence analysis or other methods.

**Test Sensitivity:**

Greater than 95% of patients with LHON are expected to have one of the variants included in this panel.<sup>1</sup> Heteroplasmy as low as 1.5% is expected to be detected by this test.

**References:**

1. Yu-Wai-Man et al. (2009) J. Med. Genet. 46 (3):145-58 (PMID: 19001017)
2. Lott et al. (2013) Curr Protoc Bioinformatics 44 :1.23.1-26 (PMID: 25489354)