

ExonArrayDx

Exon-focused array CGH analysis for detecting partial or full gene deletions or duplications

Background:

ExonArrayDx is a GeneDx-designed oligonucleotide microarray constructed for comparative genomic hybridization (oligo aCGH) to identify exonic deletions or duplications within or encompassing a targeted gene. This array targets clinically significant genes that are available on the GeneDx test menu (<http://www.genedx.com/test-catalog/genes/>), which can be ordered individually or as part of a clinically defined panel.

Custom ExonArrayDx (exome array) is a similarly designed microarray that allows targeted deletion/duplication analysis of virtually any gene(s) in the exome. Genes on the Custom ExonArrayDx can be ordered individually or as part of a custom panel of up to 20 genes.

GeneDx offers sequence analysis for a large number of metabolic and other recessive disorders. If single gene sequencing is done at GeneDx and only one mutation has been identified in a patient, and if clinically indicated, GeneDx will perform ExonArrayDx analysis to exclude a deletion of the second allele at no extra cost.

ExonArrayDx is useful as a diagnostic test in many different situations. Some common examples include:

- Testing for an autosomal dominant disorder that results from haploinsufficiency for a clinically significant gene
- Testing for an autosomal recessive disorder when sequencing or another test identifies only a single pathogenic variant
- Testing for X-linked gene deletions or duplications associated with intellectual disability

Test Methods:

ExonArrayDx probes represent the current genome build (hg19) and analysis of array data is performed using DNA Analytics software (Agilent Technologies). The array design is updated continually and optimized based on its performance and on newly published gene information. Results are compared to the HGMD database and may be confirmed by qPCR, MLPA or another array design. In cases where interpretation of results depends on whether the deletion or duplication is inherited or *de novo*, analysis of parental samples may be helpful for interpretation of the proband's results.

Test Sensitivity:

ExonArrayDx detects deletions or duplications as small as 250-500 bp. The sensitivity of ExonArrayDx is substantially higher than that of conventional testing by quantitative PCR (qPCR) or by multiplex ligation-dependent probe assay (MLPA) because it utilizes many more probes strategically positioned in or flanking most exons of a given target gene. The utility of high-resolution array CGH for single gene analysis has been demonstrated in multiple studies.¹⁻⁶

The clinical sensitivity of ExonArrayDx analysis depends on the specific gene that is being tested. Normal findings at a specific gene locus do not rule out the diagnosis of a genetic disorder associated with this gene, as a genetic abnormality may be present that is undetectable by ExonArrayDx. Specifically, Mendelian disorders that are predominantly caused by missense, nonsense, or splice site variants or by intragenic deletions or insertions that are <100bp are not detectable by ExonArrayDx. ExonArrayDx also cannot identify variants in genomic regions that are not represented on the microarray, deletions or duplications in genes that have pseudogene copies in the genome, or large unbalanced chromosomal rearrangements.

References:

1. del Gaudio D et al. Hum Mutat 29:1100-1107, 2008
2. Saillour Y et al. Hum Mutat 29:1083-1090, 2008,
3. Staaf J et al. Hum Mutat 29:555-564, 2008
4. Hegde MR et al. Hum Mutat 29:1091-1099, 2008
5. Wong LJ et al. Clin Chem 54:1141-1148, 2008
6. Dhami P et al. Am J Hum Genet 76:750-762, 2005.