

Testing for 65 Confirmed Disease-Associated Mitochondrial DNA (mtDNA) Point Variants and mtDNA Deletion Testing

Variant List: G583A, C1494T, A1555G, G1606A, G1644A, A3243G, C3256T, A3260G, T3271C, T3291C, A3302G, C3303T, G3376A, G3460A, G3635A, G3697A, G3700A, G3733A, G3733C, G3890A, C4171A, G4298A, A4300G, G4308A, G4332A, A5537insT, G5650A, G5703A, A7445G, C7471CC (=‘7472insC’), G7497A, T7511C, A8344G, T8356C, G8363A, T8528C, T8993C, T8993G, T9176C, T9176G, T9185C, T10010C, T10158C, T10191C, G10197A, T10663C, C11777A, G11778A, G12147A, G12315A, T12706C, G13051A, G13513A, A13514G, G14459A, C14482G, C14482A, T14484C, T14487C, A14495G, C14568T, T14674C, T14709C, T14849C, T14864C and large deletions

Clinical Features:

Mitochondrial disorders are clinically heterogeneous and result from dysfunction of the mitochondrial respiratory chain, which can be caused by pathogenic variants in mitochondrial DNA (mtDNA) or in nuclear genes. Some affected individuals exhibit clinical features that fall into a discrete clinical syndrome, such as Leber’s Hereditary Optic Neuropathy (LHON), chronic progressive external ophthalmoplegia (CPEO), mitochondrial encephalomyopathy with lactic acidosis and stroke-like episodes (MELAS), myoclonic epilepsy with ragged-red fibers (MERRF), neurogenic weakness with ataxia and retinitis pigmentosa (NARP) or Leigh syndrome (LS). Although the majority of mtDNA variants are rare, common mtDNA variants have been identified and are often associated with discrete clinical syndromes, whereas, other mtDNA variants have been confirmed as pathogenic because they have been described in multiple independent families. The table below lists the 65 confirmed disease-associated variants according to MITOMAP that are included in this panel.¹⁵

Genetics:

Variants in mtDNA arise de novo or are maternally inherited. In most cases, mtDNA point variants are inherited, whereas gross deletions arise de novo¹². Each mitochondrion has multiple copies of mtDNA and there are hundreds to thousands of mitochondria per cell, depending on the cell type. Usually, mtDNA variants affect only a fraction of the mtDNA; the coexistence of normal and mutant mtDNA is called heteroplasmy. When the percentage of mutant mtDNA (variant load) reaches a certain threshold that varies by tissue type, age, and specific variant the function of that tissue may become impaired.¹² As the variant load varies within and between tissues, the manifestation of mitochondrial disease may reflect the tissue-specific variant load.¹³ Many factors can affect the percent heteroplasmy these include physiologic processes that are affected by the mtDNA variant, the function of the tissue, and the rate of cell division in that tissue. Variants in mtDNA may only be identified in specific

tissues, particularly those with a lower rate of cell division such as skeletal muscle, heart and brain.¹² Large deletions of mtDNA associated with Pearson syndrome are detectable in blood, while large deletions associated with KSS and CPEO are detectable in skeletal muscle.

Test Methods:

Using genomic DNA, the entire mitochondrial genome is amplified by long-range PCR and sequenced for the detection of 65 confirmed disease-associated variants and large mtDNA deletions using a novel solid-state sequencing-by-synthesis process that allows sequencing a large number of amplicons in parallel.¹⁴ DNA sequences are assembled and compared to the published mitochondrial genome reference sequences for analysis. The presence of any variants is confirmed by conventional dideoxy sequence analysis or other methods.

Test Sensitivity:

The 65 mtDNA point variants include all mtDNA point variants that have been confirmed to be disease-associated to date. These variants account for >80% of MELAS, >80% of MERRF, >95% of LHON, >50% of MIDD, approximately 50% of NARP and approximately 20% of LS cases.¹⁵ Approximately 90% of individuals with Pearson syndrome or KSS, and 50% of patients with CPEO have a large-scale (2-10 kb) mtDNA deletion.¹⁸ MtDNA deletions larger than 2 kb account for >95% of the reported disease causing mtDNA deletions and are responsible for >99% cases of mtDNA deletion-associated mitochondrial disease.¹⁵ Overall, this test can detect pathogenic primary mtDNA defects in approximately 85% of patients. For the 65 point variants, heteroplasmy as low as 1.5% is expected to be detected and for large-scale mtDNA deletions (2 kb or larger) heteroplasmy as low as 2.5%-5% is expected to be detected. However, for large-scale mtDNA deletions observed at less than 15% heteroplasmy a quantitative value will not be provided.

Specimen Requirements and Shipping/Handling:

Special Considerations for Mitochondrial Disorders: While variants in nuclear genes are easily detectable in whole blood specimens, some mtDNA variants and deletions/duplications may only be detectable in other tissues. Tissue biopsies are preferable for mtDNA analysis, therefore, sending a blood sample together with a tissue biopsy from the same patient is recommended.

PREFERRED: TISSUE BIOPSIES (muscle or liver) AND BLOOD SPECIMEN: For tissue, please submit ≥50 mg, frozen within minutes after collection, stored at -80°C and shipped on dry ice with overnight delivery. Whole blood in EDTA; Adults: 8-10 ml; Children: 4-6 ml; Infants: 2-3 ml. Ship blood separately, overnight at ambient temperature, using a cool pack in hot weather. Blood specimens may be refrigerated for up to 7 days prior to shipping. **DO NOT FREEZE BLOOD.** **BLOOD:** Whole blood in EDTA; Adults: 8-10 ml; Children: 4-6 ml; Infants: 2-3 ml. Ship blood overnight at ambient temperature, using a cool pack in hot weather. Blood specimens may be refrigerated for up to 7 days prior to shipping. **EXTRACTED**

DNA: A minimum amount of 5 micrograms of high quality DNA, with a concentration of at least 50 ng/ul (50 nanograms per microliter). Buccal Brushes: NOT accepted for this test. Cultured fibroblasts: NOT accepted for this test

| mtDNA pathogenic variants | Examples of Associated Disorders |
|---------------------------|--|
| G583A | MELAS, Mitochondrial Myopathy and Exercise Intolerance ¹⁵ |
| C1494T | Maternally Inherited Deafness or Aminoglycoside-Induced Deafness ¹⁵ |
| A1555G | Maternally Inherited Deafness or Aminoglycoside-Induced Deafness ¹⁵ |
| G1606A | Ataxia, Myoclonus and Deafness ¹⁵ |
| G1644A | Hypertrophic Cardiomyopathy Plus MELAS ¹⁵ |
| A3243G | MELAS (3243A>G present in ~80% of cases) ¹ Maternally Inherited Diabetes and Deafness (MIDD) (3243A>G present in ~ 2%-7% of patients) ² Leigh Syndrome ¹ Hypertrophic Cardiomyopathy (3243A>G present in ~10% of Finnish patients) ² Sensorineural Hearing Loss, Focal Segmental Glomerulosclerosis, Cardiac Plus Multi-Organ Dysfunction ¹⁵ Chronic Progressive External Ophthalmoplegia / Mitochondrial myopathy ¹⁹ |
| C3256T | MELAS ¹⁵ |
| A3260G | Maternal Myopathy and Cardiomyopathy ¹⁵ |
| T3271C | MELAS (3271T>C present in ~7.5% of cases) ³ |
| T3291C | MELAS, Myopathy, Deafness plus Cognitive Impairment ¹⁵ |
| A3302G | Mitochondrial Myopathy ¹⁵ |
| C3303T | Maternal Myopathy and Cardiomyopathy ¹⁵ |
| G3376A | LHON-MELAS Overlap Syndrome ¹⁵ |
| G3460A | LHON (Together 3460G>A, 11778G>A and 14484T>C account for 95% of patients with LHON) ⁴ |
| G3635A | LHON ¹⁵ |
| G3697A | MELAS/ Leigh Syndrome/ LHON and Dystonia ¹⁵ |
| G3700A | LHON ¹⁵ |
| G3733A | LHON ¹⁵ |
| G3733C | LHON ¹⁶ |
| G3890A | Progressive Encephalomyopathy / Leigh Syndrome / Optic Atrophy ¹⁵ |
| C4171A | LHON ¹⁵ |
| G4298A | Chronic Progressive External Ophthalmoplegia/ Multiple Sclerosis ¹⁵ |
| A4300G | Maternally Inherited Hypertrophic Cardiomyopathy (MICM) ⁵ |
| G4308A | Chronic Progressive External Ophthalmoplegia ¹⁵ |
| G4332A | Encephalopathy/ MELAS ¹⁵ |
| A5537insT | Leigh Syndrome ¹⁵ |
| G5650A | Myopathy ¹⁵ |
| G5703A | Chronic Progressive External Ophthalmoplegia/ Mitochondrial Myopathy ¹⁵ |
| A7445G | Sensorineural Hearing Loss ¹⁵ |
| C7471CC (=7472insC') | Progressive Encephalopathy/ Ataxia, Myoclonus and Deafness/ Moto Neuron Disease-Like ¹⁵ |
| G7497A | Mitochondrial Myopathy/ Exercise Intolerance ¹⁵ |

| | |
|---------|--|
| T7511C | Sensorineural Hearing Loss ¹⁵ |
| A8344G | MERRF (8344A>G present in over 80% of patients) ⁶ |
| T8356C | MERRF ¹⁵ |
| G8363A | MERRF ⁶ Maternally Inherited Cardiomyopathy Plus Deafness (MICM) ^{6, 15} Autism/ Leigh Syndrome/ Ataxia Plus Lipomas ¹⁵ |
| T8528C | Infantile Cardiomyopathy ¹⁵ |
| T8993C | Leigh Syndrome (LS) (~10-20% of patients have either 8993T>C or 8993T>G) ⁷ NARP (A mutation at nucleotide 8993 is estimated to be present in 20% to greater than 50% of patients. 8993T>C is less common than 8993T>G.) ⁷ |
| T8993G | Leigh Syndrome (LS) (~10-20% of patients have either 8993T>G or 8993T>C) ⁷ NARP (A mutation at nucleotide 8993 is estimated to be present in 20% to greater than 50% of patients. 8993T>G is more common than 8993T>C.) ⁷ |
| T9176C | Leigh Syndrome (LS)/NARP (present in ~ 1-5% of patients) ⁷ / Familial Bilateral Striatal Necrosis ¹⁵ |
| T9176G | Leigh Syndrome (LS) (present in ~ 1-5% of patients) ⁷ NARP (present in ~ 1-5% of patients) ⁷ Spastic Paraplegia ¹⁵ |
| T9185C | Leigh Disease/ Ataxia Syndromes/ NARP-Like Disease ¹⁵ |
| T10010C | Progressive Encephalopathy ¹⁵ |
| T10158C | Leigh Disease ¹⁵ |
| T10191C | Leigh Disease/ Leigh-Like Disease/ Epilepsy, Strokes, Optic Atrophy and Cognitive Decline ¹⁵ |
| G10197A | Leigh Disease/ Dystonia/ Stroke/ LHON and Dystonia ¹⁵ |
| T10663C | LHON ¹⁵ |
| C11777A | Leigh Disease ¹⁵ |
| G11778A | LHON (Together 11778G>A, 3460G>A and 14484T>C account for 95% of patients with LHON. Of the three 11778G>A is the most common, present in ~70% of Caucasian patients and 90% of Asian patients) ⁴ Progressive Dystonia ¹⁵ |
| G12147A | MERFF-MELAS/ Encephalopathy ¹⁵ |
| G12315A | Chronic Progressive External Ophthalmoplegia/ Kearns Sayre Syndrome ¹⁵ |
| T12706C | Leigh Disease ¹⁵ |
| G13051A | LHON ¹⁵ |
| G13513A | MELAS (rare) ⁸ Leigh Disease/ MELAS/ LHON-MELAS Overlap Syndrome ^{15, 17} Leigh Syndrome (LS) (present in ~ 1-5% of patients) ⁷ |
| A13514G | Leigh Disease/ MELAS ¹⁵ |
| G14459A | LHON (rare) ⁹ / Leigh Disease ¹⁵ |
| C14482G | LHON ¹⁵ |
| C14482A | LHON ¹⁵ |
| T14484C | LHON (Together 14484T>C, 3460G>A and 11778G>A account for 95% of patients with LHON. ⁴ 14484T>C is the most common cause of LHON in French Canadians ¹⁰) |
| T14487C | Dystonia/ Leigh Disease/ Ataxia ¹⁵ |
| A14495G | LHON ^{15, 16} |
| C14568T | LHON ¹⁵ |

| | |
|---------|--|
| T14674C | Reversible COX Deficiency Myopathy ¹⁵ |
| T14709C | Mitochondrial Myopathy Plus Diabetes Mellitus & Deafness/ Encephalopathy ¹⁵ |
| T14849C | Exercise Intolerance / Septo-Optic Dysplasia ¹⁵ |
| T14864C | MELAS ¹⁵ |

Abbreviations:

MELAS – Mitochondrial Encephalomyopathy, Lactic Acidosis, and Stroke-like Episodes

LHON – Leber Hereditary Optic Neuropathy

MERRF – Myoclonic Epilepsy and Ragged Red Muscle Fibers

MtDNA deletion syndromes predominately consist of three overlapping phenotypes that usually occur in a single individual in a family.¹⁸ The three phenotypes are Kearns-Sayre syndrome (KSS), CPEO and Pearson Syndrome.

Characteristics of Mitochondrial DNA Deletion Syndromes

| mtDNA Deletion Syndromes | Disease Characteristics | Characteristics of mtDNA Deletions ¹⁸ |
|--------------------------|--|--|
| KSS | A triad of (1) onset < 20 y/o, (2) pigmentary retinopathy, and (3) PEO, plus at least one of the following: cardiac conduction block, cerebrospinal fluid protein concentration greater than 100 mg/dL, or cerebellar ataxia | ~90% have a large-scale 1.3-10 kb deletion usually present in all tissues, but most abundant in muscle, and often undetectable in blood cells. A deletion of 4977 bp is the most common. Over 150 deletions have been associated with KSS. Large-scale duplications have also been reported. |
| CPEO | Ptosis, ophthalmoplegia, and variably severe proximal limb weakness may be the early sign of KSS. | Deletion/duplication analysis is estimated to identify a deletion in approximately 50% of patients. Deletions are confined to skeletal muscle. |
| Pearson Syndrome | Sideroblastic anemia, exocrine pancreas dysfunction, usually fatal in infancy: children who survive the disease usually | Deletions are usually more abundant in blood than other tissue types. Deletion load gradually decreases in blood and increases in muscle as the disease evolves to PEO and KSS |

go on to develop KSS. over time.

References:

1. Longo, N. (2003) *Neurol Clin N Am* 21:817-831.
2. Majamaa et al., (1998) *Am J Hum Genet* 63:447-454.
3. Goto et al., (1991) *Biochim Biophys Acta* 1097:238-40.
4. Mackey et al., (1996) *Am J Hum Genet* 59:481-485
5. Taylor et al., (2003) *J Am Coll Cardiol* 41:1786-96.
6. DiMauro, S. *Gene Reviews* (2005) MERRF.
7. Thorburn, D. *Gene Reviews* (2006) Mitochondrial DNA-Associated Leigh Syndrome and NARP.
8. DiMauro, S. *Gene Reviews* (2005) MELAS.
9. Yu-Wai-Man, P and Chinnery, P. *Gene Reviews* (2008) Leber Hereditary Optic Neuropathy.
10. Macmillan et al., (1998) *Neurology* 50:417-22.
11. Crispim et al., (2008) *Arq Bras Endocrinol Metab* 52:1228-1235.
12. Zhu et al., (2009) *Acta Biochim Biophys Sin* 41:179-187.
13. Chinnery, P. *Gene Reviews* (2006) Mitochondrial Disorders Overview.
14. Bennett S.(2004) *Pharmacogenomics* 5:433-8.
15. MITOMAP: A Human Mitochondrial Genome Database. <http://www.mitomap.org>, 2008.
16. Achilli et al., (2012) *PLoS One* 7:e42242.
17. Pulkes et al., (1999) *Ann Neurol* 46:916-9.
18. DiMaruo, S. and Hirano, M. (Updated [May 3, 2011]). Mitochondrial DNA Deletion Syndromes In: *GeneReviews at Genetests: Medical Genetics Information Resource* (database online). Copyright, University of Washington, Seattle. 1997-2012. Available at <http://www.genetests.org>
19. Jeppesen TD,et al. (2003) *Ann Neurol* 54(1):86-92.