

## *TH* Gene Analysis in Tyrosine Hydroxylase Deficiency, Dopa-Responsive Dystonia, and Autosomal Recessive Infantile Parkinsonism

**Disorder also known as:** Autosomal Recessive Segawa Syndrome

### **Clinical Features:**

Tyrosine hydroxylase (TH) deficiency is a rare autosomal recessive movement disorder with onset typically within the first years of life. It is associated with phenotypic variability that ranges from dopa-responsive dystonia (DRD) to dopa-responsive infantile parkinsonism to infantile progressive encephalopathy that is not dopa-responsive. Additional features may include hypotonia, hypokinesia, oculogyric crises and ptosis, and autonomic signs (temperature instability, hypoglycemia). A diurnal fluctuation of symptoms may be evident. Carriers are usually asymptomatic but some have been reported with restless leg symptoms and exercise-induced stiffness. TH deficiency is typically characterized by decreased levels of homovanillic acid (HVA) and 3-methoxy-4-hydroxyphenylglycol (MHPG) with normal levels of 5-hydroxyindoleacetic acid (5-HIAA) and a decreased HVA/5-HIAA ratio in cerebrospinal fluid. There is no specific biochemical test for this disorder. The phenotype associated with TH deficient dopa-responsive dystonia may significantly overlap with DRD caused by pathogenic variants in the *GCH1* gene, but may also be more complex (DRD-plus syndrome).

### **Inheritance Pattern/Genetics:**

Autosomal Recessive

### **Test Methods:**

Using genomic DNA obtained from the submitted biological material, bi-directional sequence of the coding region, splice junctions and the cAMP response element (located at -74 to -67 bp) of the *TH* gene is analyzed. If clinically indicated, if sequencing identifies a pathogenic variant on only one allele of the *TH* gene, reflex deletion/duplication testing (ExonArrayDx) may be performed at no additional charge to evaluate for a deletion/duplication of one or more exons of this gene. Pathogenic variants found in the first person of a family to be tested are confirmed by repeat analysis using sequencing, restriction fragment analysis, or another appropriate method.

### **Test Sensitivity:**

In several small studies of patients with TH deficiency diagnosed by CSF testing, analysis of the *TH* gene identified pathogenic variants in all patients (16 patients or 32/32 *TH* alleles

characterized).<sup>2,3,5,6,7,9</sup> In one study, pathogenic variants in the *TH* gene were identified in 3/17 patients without variants in the *GCH1* gene, and in 3/7 patients with DRD-plus syndrome.<sup>8</sup>

## References:

1. Van den Heuvel et al., (1998) *Hum Genet* 102:644-6
2. Wevers et al., (1999) *J Inher Metab Dis* 22:364-373
3. Janssen et al., (2000) *Ann Hum Genet* 64:375-382
4. Furukawa et al., (2001) *Neurology* 56:260-263
5. Schiller et al., (2004) *Neurology* 63:1524-1526
6. Ribases et al., (2007) *Mol Genet Metab* 92:274-277
7. Verbeek et al., (2007) *Ann Neurol* 62:422-426
8. Clot et al., (2009) *Brain* 132: 1753-1763
9. Zafeiriou et al., (2009) *Mol Genet Metab* 97:18-20
10. Pons R et al., (2010) *Mov Disord* 25(8): 1086-90
11. Willemsen et al., (2010) *Brain* 133:1810-1822
12. Furukawa & Kish. Tyrosine Hydroxylase Deficiency. 2008 Feb 8 [Updated 2014 Jul 17]. In: Pagon RA, Adam MP, Ardinger HH, et al., editors. *GeneReviews*® [Internet]. Seattle (WA): University of Washington, Seattle; 1993-2016. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK1437/>